

BEACH CHAIR

FIXED FRAME BEACH WHEELCHAIR



The Melrose beach wheelchair has been developed to provide access to the beach and water.

Large tyres provide a low pressure footprint allowing the Melrose beach wheelchair to float over the softest sand or pebble beach.

Fabricated from aluminium with a high impact nylon mesh for the upholstery, the Melrose beach wheelchair is light enough to toss in the back of the car. The wheels can be easily removed and the back folds down so storage is no problem.

BEACH CHAIR

FIXED FRAME BEACH WHEELCHAIR

STANDARD SPECIFICATIONS

Aluminium Frame 90750

Available in 14" or 16" seat widths

Lightweight non rust aluminium frame

Stainless steel bolts and bearings

Fold down back

Flip back armrests

Adjustable foot-plate bar height

Easy release wheels and castors - split pin

Black fabric sling and tension adjustable back made from high impact, quick dry, nylon mesh

Sun canopy

5 Year Frame Warranty